

DATI TECNICI / technical sheet	SEDIA PORTANTINA 2 RUOTE
	stair chair 2 wheels

93/42/CE
2007/47/CE

MATERIALI / components

1	TELO PORTANTE <i>chair sheet</i>	POLIESTERE IMPREGNATO PVC <i>polyester spread PVC</i>
2	STRUTTURA SEDIA <i>chair frame</i>	ALLUMINO E FERRO <i>aluminium and iron</i>
3	BRACCIOLI <i>arm-rest</i>	FERRO E MANICOTTO IN PLASTICA <i>iron and plastic</i>



dimensioni aperta / open sizes
dimensioni chiusa / closed sizes
peso / weight
portata / load

910 mm(h) x 520 mm(l) 710 mm (d)
910 mm(h) x 520 mm(l) 195 mm (d)
8,5 kg
150 kg

cod : SED004